

MOLECULAR SEROLOGY PROCEDURES MANUAL

Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing

Status: Published

Document ID: 77458

DATE EFFECTIVE 01/18/2024	APPROVED BY Molecular Serology Technical Leader	PAGE 1 OF 8
------------------------------	--	----------------

Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing

1 Purpose

- 1.1 Samples are run through liquid chromatography to separate peptide markers. Marker peptides are identified by mass spectrometry.

2 Check Before Starting

- 2.1 Check the Pure Air tank. If pressure on the large gauge (right gauge) is ≤ 300 psi, **do not proceed**. Contact the Proteomics Team to change the Pure Air tank if it is less than 300 psi.



- 2.2 Check the LC/MS instrument to see how many batches are in the queue. Ensure reagents are enough for batches in queue.

- 2.3 Check all volumes in reagent bottles. If reagent bottles do not have the volumes listed below, **do not proceed**. Contact the Proteomics Team to change the reagents bottles if their volumes are less than the volumes listed below.

- 2.3.1 Methanol (50% methanol + 50% HPLC water) 200 mL
- 2.3.2 Loading Solution (phase A 2% ACN + 98% HPLC water) 200 mL
- 2.3.3 Bottle A (HPLC water) 160 mL

Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.

© NYC OFFICE OF CHIEF MEDICAL EXAMINER

Qualtrax template 040621

MOLECULAR SEROLOGY PROCEDURES MANUAL

Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing

Status: Published

Document ID: 77458

DATE EFFECTIVE 01/18/2024	APPROVED BY Molecular Serology Technical Leader	PAGE 2 OF 8
------------------------------	--	----------------

2.3.4 Bottle B (ACN) 120 mL

2.4 For **samples that need to be rerun**, follow the instructions above (2.1 to 2.3 above) and then proceed to section 3.9.

3 Liquid Chromatography – Mass Spectrometer Analysis Procedure

3.1 Retrieve the following reagents:

Cytochrome C (8 pmol/ µl) at -80°C
PCM Standard
Phase A
Acetone
Acetonitrile
Isopropanol

3.2 Prepare a working solution of Cytochrome C Master Mix sufficient to resuspend all samples (see 3.3)Work on ice or in a -20°C cold tube rack.

3.3 Input into LIMS the number of regular quant samples + the number of low quant samples divided by eight. If the number is a fraction round up to the next integer. For example, if you have five regular quant and three low quant samples you have $5 + 3/8 = 5.38$, therefore input 6 into LIMS.

Cytochrome C Master Mix Reagent Ratios (multiply volumes by the number of samples to be resuspended)	
Reagent	Volume
Phase A	93 µl
Cytochrome C (8 pmol/µl)	1.45 µl
Total volume:	94.45 µl

3.4 Resuspend samples as follows:

3.5 **REGULAR Samples:** add 94.45µl of Cytochrome C Master Mix (8 pmol/µl) to reconstitute peptides.

3.6 **LOW Concentration Samples:** add 10 µl of Cytochrome C Master Mix solution to reconstitute peptides.

3.7 Vortex.

3.8 Place in refrigerated centrifuge at 4°C and spin at 18,000 g (RCF) for 30 minutes. Record instrument and temperature in LIMS.

MOLECULAR SEROLOGY PROCEDURES MANUAL

<h3 style="margin: 0;">Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing</h3>		
Status: Published		Document ID: 77458
DATE EFFECTIVE 01/18/2024	APPROVED BY Molecular Serology Technical Leader	PAGE 3 OF 8

3.9 Pipette supernatant into LC Vial (avoid pellet if present) and store at 4°C.

3.10 Prepare the LC cleaning solution to be run alongside samples and controls:

LC Cleaning Solution	
Reagent	Volume
Acetone	2 µl
Acetonitrile	9 µl
Isopropanol	9 µl
Total Volume:	20 µl /sample

3.11 Place samples, controls, and prepared cleaning solution in Eksigent refrigerated auto sampler and fill out location of each vial in the input load plate on LIMS.

3.12 Fill in LIMS output data entry sheet and select the appropriate Acquisition Method for the sample type.

Sample Type	AcqMethod
PepCalMix, Phase A	Col1 PepCalMix Col2 PepCalMix
Cleaning Sol	Col1 Cleaning Col2 Cleaning
ENeg, Ext Pos controls, Unknowns	Col1 Unknown Col2 Unknown
Reruns (High)	Col1 Unknown Hi Col2 Unknown Hi

3.12.1 Your batch should not end on a “col 1” method. If your batch does end in a “col 1” method, add an additional phase A line at the end of the batch with a “col 2” method.

3.13 Output Plate Name should be the date (YYYYMMDD), followed by a letter identifying batch order (e.g. A, B, C, etc.) an underscore, followed by analyst initials.

3.14 Ensure the output samples are in the desired injection order before loading plate (in its default - vertical mode)

3.14.1 Positive controls should be run after unknown samples when creating batch.

3.15 Export sample batch from the output plate in LIMS.

3.16 Open LIMS created batch excel sheet in the LIMS LCMS folder and confirm it is as desired and save document as .txt format.

MOLECULAR SEROLOGY PROCEDURES MANUAL

Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing


Status: Published

Document ID: 77458


DATE EFFECTIVE 01/18/2024	APPROVED BY Molecular Serology Technical Leader	PAGE 4 OF 8
------------------------------	--	----------------



- 3.17 Open analyst software.
- 3.18 Double click Hardware Configuration → Eksigent LC and MS → Activate Profile



- 3.19 Ensure that Eksigent windows is open simultaneously (both Eksigent Control Software window and the Eksigent Autosampler window).



Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.

© NYC OFFICE OF CHIEF MEDICAL EXAMINER

Qualtrax template 040621

MOLECULAR SEROLOGY PROCEDURES MANUAL

Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing

Status: Published

Document ID: 77458

DATE EFFECTIVE 01/18/2024	APPROVED BY Molecular Serology Technical Leader	PAGE 5 OF 8
------------------------------	--	----------------

3.20 Double click Build Acquisition Batch

The screenshot shows the Analyst software interface. The left-hand menu is expanded to the 'Acquire (2)' section, where 'Build Acquisition Batch' is circled in red. The main window displays a table of acquisition data with columns for Start Time, Sample Name, and Plate. A progress bar at the top indicates 100% completion for 'Acquiring Sample 20 of 20'.

	Start Time	Sample Name	Pla
1	6/14/2019 3:39:40 P	Phase A_Col2	1
2	6/14/2019 4:12:40 P	CytoC Dig_10fmol/ul 10ul>Loading_Col1	1
3	6/14/2019 4:45:41 P	CytoC Dig_10fmol/ul 10ul>Loading_Col2	1
4	6/14/2019 5:18:40 P	CytoC Dig_10fmol/ul 10ul>Loading_Col1	1
5	6/14/2019 5:51:41 P	Phase A_cbl2	1
6	6/14/2019 6:24:40 P	Phase A C1	1
7	6/14/2019 7:14:24 P	Phase A C2	1
8	6/14/2019 8:04:06 P	SE 0.1ug/ul, 3 ul loading, Col1_Sample 1	1
9	6/14/2019 8:53:49 P	SE 0.1ug/ul, 3 ul loading, Col1_Sample 1	1
10	6/14/2019 9:43:33 P	SE 0.1ug/ul, 3 ul loading, Col1_Sample 1	1
11	6/14/2019 10:33:14	SE 0.1ug/ul, 3 ul loading, Col1_Sample 1	1
12	6/14/2019 11:22:55	SE 0.1ug/ul, 3 ul loading, Col1_Sample10	1
13	6/15/2019 12:12:38	SE 0.1ug/ul, 3 ul loading, Col1_Sample 1	1
14	6/15/2019 1:02:20 A	SE 0.1ug/ul, 3 ul loading, Col1_Sample 1	1

3.21 Right click to import acquisition batch.

The screenshot shows the Analyst software interface. The left-hand menu is expanded to the 'Acquire (7)' section, where 'Build Acquisition Batch' is circled in red. The main window displays a 'Batch Editor' dialog box with a table of acquisition data. A context menu is open over the table, with 'Import From' circled in red.

Sample Name	Batch Code	Batch Position	Plate Code	Plate Position	Vial Position
Open...					
Import From					
Save As Batch					
Save As a Template...					
Hide/Show Column...					
Save Column Settings...					

Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.

© NYC OFFICE OF CHIEF MEDICAL EXAMINER

Qualtrax template 040621

MOLECULAR SEROLOGY PROCEDURES MANUAL


Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing

Status: Published


Document ID: 77458

DATE EFFECTIVE 01/18/2024	APPROVED BY Molecular Serology Technical Leader	PAGE 6 OF 8
------------------------------	--	----------------

3.22 Select .txt batch list exported from LIMS. Click open.



3.23 Select autosampler (ekspert nanoLC 400). Click OK.



3.24 **Label and Position WITNESS:** Have a witness verify the selected autosampler, batch sample names, methods, and tube positions in autosampler match that in Vial Position column.

Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.

© NYC OFFICE OF CHIEF MEDICAL EXAMINER

Qualtrax template 040621

MOLECULAR SEROLOGY PROCEDURES MANUAL


Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing

Status: Published


Document ID: 77458

DATE EFFECTIVE 01/18/2024	APPROVED BY Molecular Serology Technical Leader	PAGE 7 OF 8
------------------------------	--	----------------


3.25 Click on Submit tab when sample list is ready.



3.26 Click Submit button.



3.27 Ensure "All samples" is selected and that the box for "apply new samples order" is unchecked.



Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.


© NYC OFFICE OF CHIEF MEDICAL EXAMINER

Qualtrax template 040621


MOLECULAR SEROLOGY PROCEDURES MANUAL

<h3 style="margin: 0;">Body Fluid Identification by Proteomic Mass Spectrometry - Liquid Chromatography & Mass Spectrometer Processing</h3>		Document ID: 77458
Status: Published		
DATE EFFECTIVE 01/18/2024	APPROVED BY Molecular Serology Technical Leader	PAGE 8 OF 8


3.28 Click Queue button and double check all samples were submitted.



3.29 Click Ready button.



3.30 Click Start sample button. The Eksigent LC and 6500 MS will process all samples.



Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.

© NYC OFFICE OF CHIEF MEDICAL EXAMINER

Qualtrax template 040621