

Amplification using the Power Plex Fusion System		
Status:Published		Document ID: 5980
DATE EFFECTIVE 02/02/2018	APPROVED BY Nuclear DNA Technical Leader	PAGE 1 OF 4

PowerPlex® Fusion Sample Preparation for Amplification

1 Procedure

Fusion Sample Input Amount
Optimal – 525pg*
Minimum – 37.5pg

*The option for amplification with a greater input amount is available if determined appropriate for the sample by the analyst.

- 1.1 Retrieve the following reagents from the associated refrigerator and/or freezer.

PowerPlex Fusion® 5X Primer Pair Mix
PowerPlex Fusion® 5X Master Mix
Water, Amplification Grade
2800M Control DNA, .250ng/µl

- 1.2 Retrieve sample(s) needed for amplification from associated refrigerator and/or freezer.
- 1.3 Prepare dilutions **using Promega Amplification Grade Water**, for each sample, if necessary, according to Table 1. Vortex and centrifuge samples prior to aliquoting for dilution.

TABLE 1: Dilutions

Dilution	Amount of DNA Template (uL)	Amount of Promega® Water (uL)
0.25	3 or (2)	9 or (6)
0.2	2	8
0.1	2	18
0.05	2	38
0.04	4 or (2)	96 or (48)
0.02	2 or (1)	98 or (49)
0.01	2	198
0.008	4 or (2)	496 or (248)

- 1.4 Centrifuge reagent tubes briefly to bring contents to the bottom and then vortex for 15 seconds before use. Do NOT re-centrifuge the Master Mix or Primer Pair Mix as this may cause the reagents to be concentrated at the bottom of the tube.

FORENSIC BIOLOGY PROTOCOLS FOR FORENSIC STR ANALYSIS

Amplification using the Power Plex Fusion System		
Status:Published		Document ID: 5980
DATE EFFECTIVE 02/02/2018	APPROVED BY Nuclear DNA Technical Leader	PAGE 2 OF 4

- 1.5 Consult the Reagents tab in LIMS for the exact amount of PowerPlex Fusion® 5X Primer Pair Mix and PowerPlex Fusion® 5X Master Mix to add.

Reagent	Per reaction
5X Primer Pair Mix	2.5 µL
5X Master Mix	2.5 µL
Mastermix total:	5 µL
DNA	7.5 µL

- 1.6 Vortex prepared Master Mix and samples created in Step 3. After vortexing, **briefly centrifuge** master mix and samples.
- 1.7 Add **5 µL** of the prepared master mix to each tube that will be utilized, changing pipette tips and remixing master mix as needed.
- 1.8 **Witness Step.** Have another analyst witness the sample set-up.
- 1.9 Positive Control – total input amount of 500pg.
- 1.9.1 Aliquot positive control according to amplification sheet
- 1.10 Amplification Negative
- 1.10.1 7.5 uL of Water, Amplification Grade
- 1.11 Samples
- 1.11.1 Aliquot samples according to amplification sheet
- 1.12 Ensure that all caps are properly closed prior to sending the samples to the post-amplification laboratory.
- 1.13 Spin down samples at 1000 RPM for one minute.

2 PowerPlex® FusionPCR Conditions for the Applied Biosystems GeneAmp PCR System 9700

- 2.1 Turn on the ABI 9700 Thermal Cycler.
- 2.2 Choose the following program in order to amplify these samples:

PowerPlex® Fusion
user: casework
file: PPFusion-29

Controlled versions of Department of Forensic Biology Manuals only exist in the Forensic Biology Qualtrax software. All printed versions are non-controlled copies.

FORENSIC BIOLOGY PROTOCOLS FOR FORENSIC STR ANALYSIS

Amplification using the Power Plex Fusion System		
Status: Published		Document ID: 5980
DATE EFFECTIVE 02/02/2018	APPROVED BY Nuclear DNA Technical Leader	PAGE 4 OF 4

2.9.2 Fill out the performed by tab for the Amplification Run Review.

2.9.3 Select all output samples and click Review to perform the test batch approval.

2.10 Schedule the samples to the appropriate STR test batch and create the test batch.

ARCHIVED