

FORENSIC BIOLOGY PROTOCOLS FOR FORENSIC STR ANALYSIS

Reconstitution of Evaporated Extracts		
Status: Published		Document ID: 48256
DATE EFFECTIVE 06/30/2021	APPROVED BY Nuclear DNA Technical Leader	PAGE 1 OF 1

Reconstitution of Evaporated Extracts

1 LIMS Pre-Processing

- 1.1 Samples not already present in LIMS will need to be added into the system. See [Creating Ad Hoc Samples](#).

2 Procedure

- 2.1 Retrieve samples.
- 2.2 Retrieve and label one 1.5mL tube as a reconstitution negative.
- 2.3 **Witness Step:** Have a witness confirm the physical tube top and sample labels and the sample order against the listed samples in LIMS.
- 2.4 Using a pipette, measure the volume (if any) of liquid present in the evaporated sample(s). Log the volume.
- 2.5 Add an amount of 0.1X TE buffer necessary to bring the total volume of each sample to the required volume. Log the volume added for each sample.
- 2.6 Vortex and centrifuge each sample briefly.
- 2.7 Submit samples for quantification.

3 LIMS Post-Processing

- 3.1 Complete the data entry for all samples within LIMS.
- 3.2 In the data entry screen, all samples can be assigned a storage location in the next available storage box.
- 3.3 Store the extracts at 2° to 8°C, or frozen.